Produktinformation
Product Information
Données du produit
Instrucciones de uso
Informazioni sul prodotto

IMPRIMO® LC Model

IMPRIMO® LC Splint

IMPRIMO® LC Impression

IMPRIMO® LC Cast

IMPRIMO® LC Denture

IMPRIMO® LC Gingiva

IMPRIMO® LC Temp

IMPRIMO® LC Temp It

IMPRIMO® LC IBT

IMPRIMO® LC Try-In

IMPRIMO® LC MJF

IMPRIMO® LC Model Produktinformation

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C € 0044

IMPRIMO® LC Model Product Information

Description:

IMPRIMO® LC Model is a light-curing resin based on methacrylate and specifically engineered for the fabrication of dental models; photo-polymerization is done by means of image projection systems (≤ 405 nm). The printed models are characterized by similar haptical and optic qualities like conventional plaster models.

IMPRIMO® LC Model should be used in building processes with a z-resolution of 10-100 μm .

Processing instructions/application:

Make sure IMPRIMO® LC Model is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity.

Leave the material to stand for approx. one hour to avoid any air inclusions.

Contra-indications:

IMPRIMO® LC Model should not be used for any purpose other than the generative CAD / CAM fabrication of dental models and is not suitable for application in the mouth.

IMPRIMO® LC Model is a laboratory product.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required.

The polymerization process in the light oven in a nitrogen environment takes 5 minutes. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Model might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC Model. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.**Technical data:**

Density:	approx. 1.1 g/m
Viscosity:	approx. 0.7 Pa s
Cured material:	(depends on irradiation unit)
Flexural strength:	108 Mpa
Elongation at break:	5 %
Shore hardness:	85 D

* These data result from measurements of a representative sample.

Contains:

alkoxylated bisphenol-A dimethacrylate, initiators, stabilizers, dyes and pigments.

IMPRIMO® LC Splint Product Information

Description:

IMPRIMO® LC Splint is a light-curing resin based on methacrylate and specifically engineered for the fabrication of high-precision occlusal splints and surgical guides; photo-polymerization is done by means of image projection systems (≤ 405 nm).

IMPRIMO® LC Splint has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC Splint should be used in building processes with a z-resolution of 25-100 μm . A minimum wall thickness of 1.5 mm is required for reasons of stability.

Processing instructions/application:

Make sure IMPRIMO® LC Splint is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions.

Contra-indications:

IMPRIMO® LC Splint should not be used for any purpose other than the generative CAD/CAM fabrication of occlusal splints and surgical guides. IMPRIMO® LC Splint is not suitable for printing aligners. IMPRIMO® LC Splint is a medical device class IIa.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required. The polymerization process in the light oven takes 3 minutes in a nitrogen environment. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO LC Splint might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC Splint. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin

contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.

Technical data:

Density: approx. 1.1 g/m

Viscosity: approx. 0.7 Pa s

Cured material: (depends on irradiation unit)

Flexural strength: 93 MPa

Elongation at break: 10 %

Solubility: 1.4 $\mu\text{g mm}^{-3}$

Water absorption: 24 $\mu\text{g mm}^{-3}$

Shore hardness: 80 D

Comply: biocompatibility (irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity)

* These data result from measurements of a representative sample.

Contains:

alkoxylated bisphenol-A dimethacrylate, initiators, stabilizers, dyes and pigments.

CE 0044

IMPRIMO® LC Impression Product Information

Description:

IMPRIMO® LC Impression is a light-curing resin based on methacrylate and specifically engineered for the fabrication of high-precision individual trays; photo-polymerization is done by means of image projection systems (≤ 405 nm).

IMPRIMO® LC Impression has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC Impression should be used in building processes with a z-resolution of 50-300 μ m.

Processing instructions/application:

Make sure IMPRIMO® LC Impression is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions.

Contra-indications:

IMPRIMO® LC Impression should not be used for any purpose other than the generative CAD/CAM fabrication of individual impression trays.

IMPRIMO® LC Impression is a medical product class I.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required. The polymerization process in the light oven in a nitrogen environment takes 5 minutes. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Impression might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC Impression. In the event of eye contact, rinse thoroughly with water and seek medical

advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.

Technical data:

Density: approx. 1.1 g/m

Viscosity: approx. 0.7 Pa s

Cured material: (depends on irradiation unit)

Flexural strength: 84 MPa

Elongation at break: 10 %

Shore hardness: 80 D

Comply: biocompatibility (irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity)

* These data result from measurements of a representative sample.

Contains:

alkoxylated bisphenol-A dimethacrylate, initiators, stabilizers and dyes.

IMPRIMO® LC Cast Product Information

Description:

IMPRIMO® LC Cast is a light-curing resin based on methacrylate and specifically engineered for the distortion-free fabrication of castings; photo-polymerization is done by means of image projection systems (≤ 405 nm). Processing with standard investment materials.

IMPRIMO® LC Cast should be used in building processes with a z-resolution of 10-100 μm .

Processing instructions/application:

Make sure IMPRIMO® LC Cast is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity.

Leave the material to stand for approx. one hour to avoid any air inclusions.

Contra-indications:

IMPRIMO® LC Cast should not be used for any purpose other than the generative CAD / CAM fabrication of castings.

IMPRIMO® LC Cast is a laboratory product.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required. The polymerization process in the light oven in a nitrogen environment takes 5 minutes. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit.

Any deviation from the described manufacturing process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Cast might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC Cast. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.
Technical data:

Density:	approx. 1.1 g/m
Viscosity:	approx. 0.3 Pa s
Cured material:	(depends on irradiation unit)
Flexural strength:	86 MPa
Elongation at break:	5 %
Shore hardness:	85 D

* These data result from measurements of a representative sample.

Contains:

urethane dimethacrylates, initiators, stabilizers and dyes.

IMPRIMO® LC Denture Product Information

Description:

IMPRIMO® LC Denture is a light-curing resin based on methacrylate and specifically engineered for the fabrication of denture bases; photo-polymerization is done by means of image projection systems (≤ 405 nm).

IMPRIMO® LC Denture has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC Denture should be used in building processes with a z-resolution of 50- 100 μ m.

Processing instructions/application:

Make sure IMPRIMO® LC Denture is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions. Intraoral contact with uncured or partially uncured material shall be excluded.

The following specifications shall be considered for the digital design of the print object:

Alveolar ridge lingual ≥ 2.5 mm (upper / lower jaw)

Palatal / lingual ≥ 3.0 mm (upper jaw) ≥ 2.5 mm (lower jaw)

Vestibular / buccal ≥ 2.5 mm (upper / lower jaw)

We recommend using the bonding agent VITA VIONIC® BOND for fixing prosthetic teeth in the denture base (for ex. teeth printed using IMPRIMO® LC Resin). It may be required to additionally roughen the surfaces. The generally accepted professional techniques and dental procedures shall be carefully followed when preparing dentures.

Contra-indications:

IMPRIMO® LC Denture should not be used for any purpose other than the generative CAD / CAM fabrication of denture bases.

IMPRIMO® LC Denture is a medical product class IIa.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required. The polymerization process in the light oven in a nitrogen environment takes 10 minutes. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Any deviation from the manufacturing method described or storage conditions may affect the mechanical and optical material properties.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Denture might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC Denture. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

May cause long-lasting harmful effects to aquatic life.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.

Technical data:

Viscosity: approx. 0.5 Pa s

Cured material: (depends on irradiation unit)

Flexural strength: 114 MPa

Elongation at break: 8 %

Solubility: 0.41 µg mm⁻³

Water absorption: 25.8 µg mm⁻³

Shore hardness: 85 D

Comply: biocompatibility (irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity)

* These data result from measurements of a representative sample.

Contains:

functional methacrylate resins, initiators, dyes and stabilizers.

CE 0044

IMPRIMO® LC Gingiva Product Information

Description:

IMPRIMO® LC Gingiva is a light-curing resin based on methacrylate and specifically engineered for the fabrication of gingival masks; photo-polymerization is done by means of image projection systems (≤ 405 nm).

IMPRIMO® LC Gingiva should be used in building processes with a z-resolution of 25-100 μm .

Processing instructions/application:

Make sure IMPRIMO® LC Gingiva is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Dental molded parts can be coated and repaired in the conventional way.

Leave the material to stand for approx. one hour to avoid any air inclusions.

Contra-indications:

IMPRIMO® LC Gingiva should not be used for any purpose other than the generative CAD/CAM fabrication of gingival masks.

IMPRIMO® LC Gingiva is a laboratory product.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required. The polymerization process in the light oven in a nitrogen environment takes 5 minutes. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Gingiva might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC Gingiva. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.**Technical data:**

Density:	approx. 1.1 g/m
Viscosity:	approx. 2.5 Pa s
Cured material:	(depends on irradiation unit)
Elongation at break:	90%
Tear strength:	5 MPa
Shore hardness:	50 A

* These data result from measurements of a representative sample.

Contains:

urethane methacrylates, initiators, stabilizers and dyes.

IMPRIMO® LC Temp Product Information

Description:

IMPRIMO® LC Temp is a light-curing resin based on methacrylate and specifically engineered for the fabrication of temporary frontal restorations; photo-polymerization is done by means of image projection systems (≤ 405 nm).

IMPRIMO® LC Temp has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC Temp should be used in building processes with a z-resolution of 50-100 μm .

Processing instructions/application:

Make sure IMPRIMO® LC Temp is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity.

Leave the material to stand for approx. one hour to avoid any air inclusions.

Temporary restorations can be individualized using light-curing composites. The reduction of the restoration should not exceed 0.3 mm (incisal and vestibular areas). Please follow the instructions of the manufacturer of light-cured composite materials. For cementation, common temporary tooth cements (that do not contain eugenol) and adhesives can be used. The inner areas of the crown should be conditioned by sandblasting or other commonly used techniques.

Contra-indications:

IMPRIMO® LC Temp should not be used for any purpose other than the generative CAD/CAM fabrication of temporary frontal restorations. Bridges with a maximum of one pontic can be fabricated. When realising frontal bridges, ensure that the connection regions must have at least 12 mm² and in the case of posterior bridges at least 14 mm². The occlusal wall thicknesses must be at least 1.5 mm (central fissure) and the circumferential thickness 1 mm.

IMPRIMO® LC Temp is a medical product class IIa.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required. The polymerization process in the light oven in a nitrogen environment takes 10 minutes. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit.

Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Temp might occur, for more details please refer to our safety data sheet on IMPRIMO® LC Temp. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

May cause long-lasting harmful effects to aquatic life.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.

Technical data:

Viscosity: approx. 1.5 Pa s

Cured material: (depends on irradiation unit)

Flexural strength: 113 MPa

Elongation at break: 4 %

Solubility: 1.1 µg mm⁻³

Water absorption: 31.3 µg mm⁻³

Shore hardness: 80 D

Comply: biocompatibility (irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity)

* These data result from measurements of a representative sample.

Contains:

functional methacrylate resins and inorganic fillers with particle sizes from 0.4 to 3 µm.

CE 0044

IMPRIMO® LC Temp It Product Information

Description:

IMPRIMO® LC Temp It is a light-curing resin based on methacrylate and specifically engineered for the fabrication of long-term temporaries such as crowns, bridges and mock-ups; photo-polymerization is done by means of image projection systems (≤ 405 nm).

IMPRIMO® LC Temp It has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC Temp It should be used in building processes with a z-resolution of 50-100 μm .

Processing instructions/application:

Make sure IMPRIMO® LC Temp It is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity.

Leave the material to stand for approx. one hour to avoid any air inclusions.

Temporary restorations can be individualized using light-curing composites. The reduction of the restoration should not exceed 0.3 mm (incisal and vestibular areas). Please follow the instructions of the manufacturer of light-cured composite materials. For cementation, common temporary tooth cements (that do not contain eugenol) and adhesives can be used. The inner areas of the crown should be conditioned by sandblasting or other commonly used techniques.

Contra-indications:

IMPRIMO® LC Temp It should not be used for any purpose other than the generative CAD / CAM fabrication of long-term temporaries such as crowns, bridges and mock-ups. Bridges with a maximum of one pontic can be fabricated. When realising frontal bridges, ensure that the connection regions must have at least 12 mm² and in the case of posterior bridges at least 14 mm². The occlusal wall thicknesses must be at least 1.5 mm (central fissure) and the circumferential thickness 1 mm.

IMPRIMO® LC Temp It is a medical product class II a.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required. The polymerization process in the light oven in a nitrogen environment takes 10 minutes. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Temp It might occur, for more details please refer to our safety data sheet on IMPRIMO® LC Temp It. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

May cause long-lasting harmful effects to aquatic life.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.

Technical data:

Viscosity: approx. 3.5 Pa s

Cured material: (depends on irradiation unit)

Flexural strength: 169 MPa

Elongation at break: 4 %

Solubility: 1.1 µg mm⁻³

Water absorption: 31.1 µg mm⁻³

Shore hardness: 80 D

Comply: biocompatibility (irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity)

* These data result from measurements of a representative sample.

Contains:

functional methacrylate resins and inorganic fillers with particle sizes from 0.4 to 3 µm.

CE 0044

IMPRIMO® LC IBT Product Information

Description:

IMPRIMO® LC IBT is a light-curing resin based on methacrylate and specifically engineered for the fabrication of a transfer matrix in the indirect bonding technique; photo-polymerization is done by means of image projection systems (≤ 405 nm).

IMPRIMO® LC IBT has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC IBT should be used in building processes with a z-resolution of 50-100 μm .

Processing instructions/application:

Make sure IMPRIMO® LC IBT is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions.

Contra-indications:

IMPRIMO® LC IBT should not be used for any purpose other than the generative CAD/CAM fabrication of transfer matrixes in the indirect bonding technique.

IMPRIMO® LC IBT is a medical product class I.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required. The polymerization process in the light oven in a nitrogen environment takes 7 minutes. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC IBT might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC IBT. In the event of eye

contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.

Technical data:

Density: 1.1g/ml

Viscosity: approx. 2.5 Pa s

Cured material: (depends on irradiation unit)

Elongation at break: 50%

Tear strength: 6, 2MPa

Shore hardness: 40 D

Comply: Biocompatibility (irritation and delayed-type hypersensitivity; cytotoxicity)

* These data result from measurements of a representative sample.

Contains:

urethane methacrylates, initiators, stabilizers and dyes.

IMPRIMO® LC Try-In Product Information

Description:

IMPRIMO® LC-Try-in is a light-curing resin based on methacrylate and specifically engineered for the fabrication of individual try-ins; photo-polymerization is done by means of image projection systems (≤ 405 nm).

IMPRIMO® LC-Try-In has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC- Try-In should be used in building processes with a z-resolution of 50-100 μm .

Processing instructions/application:

Make sure IMPRIMO® LC-Try-In is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity.

Dental molded parts can be repaired conventionally.

Leave the material to stand for approx. one hour to avoid any air inclusions.

Contra-indications:

IMPRIMO® LC-Try-In should not be used for any purpose other than the generative CAD / CAM fabrication of individual try-ins.

IMPRIMO® LC-Try-In is a medical product class I.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required. The polymerization process in the light oven in a nitrogen environment takes 5 minutes. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Try-In might occur, for more details please refer to our safety data sheet on IMPRIMO® LC Try-In.

In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.

Technical data:

Density: approx. 1.1 g/m

Viscosity: approx. 0.7 Pa s

Cured material: (depends on irradiation unit)

Flexural strength: 91 MPa

Elongation at break: 8 %

Shore hardness: 80 D

Comply: biocompatibility (irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity)

* These data result from measurements of a representative sample.

Contains:

alkoxylated bisphenol-A dimethacrylate, initiators, stabilizers and dyes.

IMPRIMO® LC MJF Product Information

Description:

IMPRIMO® LC MJF is a light-curing resin based on methacrylate and specifically engineered for the fabrication of a devices for oral and maxillo-facial surgery as well as for implantology; photo-polymerization is done by means of image projection systems (≤ 405 nm). Suitable for current sterilization methods such as plasma sterilization, gamma sterilization, ethylene oxide sterilization and autoclave sterilization.

IMPRIMO® LC MJF has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC MJF should be used in building processes with a z-resolution of 25-100 μ m.

Processing instructions/application:

Make sure IMPRIMO® LC MJF is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions.

Contra-indications:

IMPRIMO® LC MJF should not be used for any purpose other than the generative CAD / CAM fabrication of devices for oral and maxillo-facial surgery as well as implantology.

IMPRIMO® LC MJF is a medical product class IIa.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post-exposure:

Final and complete curing of the objects is done in the curing unit IMPRIMO® Cure. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of ≤ 405 nm is required. The polymerization process in the light oven in a nitrogen environment takes 5 minutes. The corresponding parameters are stored in the internal programme of the IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC MJF

might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC MJF. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention:

Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Printer parameters available for Asiga printers.

For professional use only.

Technical data:

Density:	1.1 g/ml
Viscosity:	approx. 0.7 Pa s
Cured material:	(depends on irradiation unit)
Flexural strength:	117 MPa
Elongation at break:	5 %
Solubility:	0.51 $\mu\text{g mm}^{-3}$
Water absorption:	19.9 $\mu\text{g mm}^{-3}$
Shore hardness:	85 D
Comply:	biocompatibility (hemocompatibility; irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity).

* These data result from measurements of a representative sample.

Contains:

alkoxylated bisphenol-A dimethacrylate, initiators, stabilizers and dyes.

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SCHEU-DENTAL GmbH